




Use it anywhere



No water is used because the waste is sealed using heat compression. It can be installed anywhere by connecting it with a car DC cable.

 ${\rm *The~unit~comes~standard~with~an~AC~power~supply.~*Can~be~connected~using~the~optional~car~DC~cable~CA-50.}$



Lightweight and compact

The total weight is approximately 7.0 kg, making it lightweight and easy to move. It is small and easy to store.

It doesn't take up much space and comes with a storage bag for easy transport.





Easy post-processing



The film used in these wrap bags is an environmentally friendly biodegradable film. The sealed wrap bag can be disposed of in the same way as a disposable diaper.







USE Easy to use



Add special coagulant

A special coagulant is administered and the excretion is completed. A coagulant is used to solidify the water in the excrement. Use.



Press the switch button

After excretion, press the switch button. The lap will start automatically. The film is fed automatically.



Remove the wrapped bag

After about 100 counts, the end sound (beep beep beep) When it beeps, the disconnection is complete. Remove the wrapped bag.

Wrappon Sunny

Dimensions	Closed D46.9 x W36 x H42.1cm
	folded D46.5 x W36 x H24.7cm
Weight	7.0kg
Power	100-240V, 50/60Hz, 60w jmax consumption
Materials	Stainless steel, ABS,PP,POM
Load capacity	150kg
Included items	Consumables : 30uses(film rolls and coagulant) x 2sets. Car DC Calbe



Contact information

